METTLER TOLEDO

Smart sensor support you to define quality

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Weight serves a critical process parameter or quality attribute in many Up- and Downstream unit operations

Up-stream Processing - Fermentation

- Control and measurement of feeds like media, buffer, antifoaming agents
- Measure reactor content

Downstream Processing

- Automated dosing of reagents
- Accurate Formulations

In Media/Buffer Preparation

- Precise Pre-weighing of ingredients
- Mixing and blending in tanks

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
Automation Transactions for Error Free Production

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(41) BALANCES USP34

Change to read:

■ This chapter states the requirements for balances used for materials that must be accurately weighed (see *General Notices*, 8.20). Unless otherwise specified, when substances must be "accurately weighed", the weighing shall be performed using a balance that is calibrated over the operating range and meets the requirements defined for repeatability and accuracy. For balances used for other applications, the balance repeatability and accuracy should be commensurate with the requirements for its use.

For discussion of the theoretical basis of these requirements, see general information chapter *Weighing on an Analytical Balance* (1251), which may be a helpful—but not mandatory—resource.

REPEATABILITY

Repeatability is assessed by weighing one test weight NLT 10 times. [NOTE—The test weight must be within the balance's operating range, but the weight need not be cali-

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EP

CALIBRATION

Calibration is part of balance qualification and is performed by the user or by a suitable competent body. Its aim is to establish traceability of measurement results to SI units (metrological traceability). The calibration results include measurement uncertainty and are documented in a calibration certificate. To ensure traceability, it is recommended to perform calibration before any maintenance operation is carried out on the balance that significantly alters its measurement performance. 'Significant' operations include repairs, transfer of the balance to another location or mechanical adjustment of one or more weighing parameters. The balance must be re-calibrated after significant operations. Re-calibration is not necessary after less significant operations, which include levelling the balance or adjustments using built-in weights.

PERFORMANCE CHECKS

Performance checks are carried out to evaluate the random and systematic error of a balance; they consist of measuring precision and accuracy respectively and comparing the results obtained to pre-defined acceptance criteria. Balances are considered suitable if none of these errors exceeds 0.10 per cent.


In practice, performance checks focus on the two weighing parameters that most significantly affect the performance of the instrument, i.e. repeatability, for precision, and sensitivity as the main component of the accuracy of the balance.

Accuracy is also impacted by two other parameters: eccentricity and linearity. A quadratic addition of the errors of these individual parameters, rather than a more conservative linear addition, provides a more realistic approach to the assessment of the accuracy of the balance because the three individual parameters are known to be largely independent of each other, and it is considered unlikely that they will occur simultaneously and have the same algebraic sign. Therefore, the acceptance criterion for each individual parameter can be set at 0.05 per cent, i.e. half the overall accuracy tolerance of 0.10 per cent. While accuracy is impacted by all three parameters, the impact of eccentricity and linearity is typically less than that of sensitivity. Hence,

"Operation that, under specified conditions, in a first step, establishes a relation between the quantity values **with measurement uncertainties** provided by measurement standards and corresponding indications with **associated measurement uncertainties** and, in a second step, uses this information to establish a relation for obtaining a measurement result from an indication."

International Vocabulary Of Basic And General Terms In Metrology (VIM)
JCGM 200:2012, item 2.39

In simple words: During calibration, the measurement uncertainty of the instrument is determined.

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
Calibration Certificate ID
TMS643-04-022223-ACC

METTLER TOLEDO

Service Provider
3F No. 17, Lane 171, Ju-Zong Rd., Sec.2 Taipei 1149
NA
NA

Customer
Accuracy Calibration Certificate

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Weighing Device

Manufacturer: Mettler Toledo	Instrument Type: Weighing Instrument
Model: PB3604PW-CC300	Asset Number: S1 FEAJ SGS
Serial No.: C011152843	Transfer Model: N/A
Year: 1978	Transfer Model No.: N/A
Plant: ZF	Transfer Asset No.: N/A
Room: N/A	Transfer Asset No.: N/A

Range	Max. Capacity	Resolution (g)
1	300 kg	0.0005 kg

Procedure

Calibration Standard: EURAMET (cp-18) = 4.0 (11/2018)
METTLER TOLEDO Work Instruction: General_SOP_ACC_30000553

This calibration certificate contains measurements for As Found calibration. No As Left calibration was performed because the device was not modified after As Found calibration. Therefore, results for As Left correspond to As Found.

In accordance with EURAMET (cp-18) (11/2018), the test loads were selected to reflect the specific use of the weighing device or to accommodate specific calibration conditions.

As Found Calibration Date: 10-02-2023
As Left Calibration Date: N/A
Issue Date: 22-02-2023
Next Calibration Date: 10-02-2024

Service Technician: Stanley Lin
Stanley Lin

Software Version: 7.02.1.79
Report Version: 1.0.0.0
Print Number: 002_000

© METTLER TOLEDO
This is an original document and may not be printed, reproduced or altered without the written permission of the issuing calibration laboratory.

Calibration Certificate ID
TMS643-04-022223-ACC

METTLER TOLEDO Service

Measurement Uncertainty of the Weighing Instrument in Use

Stated in the expanded uncertainty with $k=2$ in use. The formula shall be used for the estimation of the uncertainty under consideration of the errors of indication. The value R represents the net load indication in the unit of measure of the device.

Temperature coefficient for the evaluation of the measurement uncertainty in use: $3.0 \cdot 10^{-4} / K$
Temperature range in use for the evaluation of the measurement uncertainty in use: $4 K$

Uncertainty of Uncertainty Budget

Range	U	As Found	As Left
1 0.0005 kg - 300 kg	$U_1 = 0.86 g + 0.0708 g/kg \cdot R$		N/A

To optimize the stability of the instrument, besides of the zero load only, increasing measurement points with a net load of 5% of measurement range or larger are suggested. The user equipment shall be used in the same way as during the calibration.

The linear approximation of the uncertainty in use should be used as a breakdown only. The underlying distribution of the errors of indication does not stipulate a static approximation function (e.g., μ and σ and equation).

Absolute and Relative Measurement Uncertainty in Use for Various Net Load Values (Examples)

Net Indication	As Found	As Found	As Left	As Left
0.0000 kg	0.86 g	0.1%	N/A	N/A
0.0000 kg	0.86 g	0.1%	N/A	N/A
3.0000 kg	1.1 g	0.037%	N/A	N/A
30.0000 kg	3.0 g	0.010%	N/A	N/A
300.0000 kg	23 g	0.0077%	N/A	N/A

Measurement Uncertainty of the Weighing Instrument in Use

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
GWP® Recommendation

3Q Pacs

StarterPac Installation
StarterPac
The StarterPac is the entry solution for all other non- or low regulated industries
IPac Standard Qualification
IPac
The IPac is for industries with high quality standards such as Pharma, Chemicals, Food & Beverage, Automotive etc.
EQPac Comprehensive Qualification
EQPac
The EQPac provides detailed lifecycle qualification documentation and is for highly regulated industries such as Pharma, Biotech, Cosmetics and Food

ACC (Uncertainty)

METTLER TOLEDO Accuracy Calibration Certificate
GWP® Certificate
Process Requirements
Scale Weighing Range

GWP® Verification

Good Weighing Practice™ GWP® Verification
Table with columns: Parameter, Value, Unit, Note

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PBK9

The Monobloc load cell at the core of the PBK9 weighing platform guarantees the highest precision and reliability. The robust load-cell housing features integrated overload protection and durable mechanical interfaces for many years of intensive use.

PB/FK



Figure 1 – View into a Weigh Module with POWERCELL®

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

Small Volume Filling and Dispensing

1. Accurate sensor confirms correct filling

- Feedback to filling device
- Flexibility in cell designs
- No contact with electrolyte
- Consistent cell quality

2. Weight-based filling

- Highest precision with best measurement technology
- Constant adjustment of filler
- Self-calibration and adjustment
- Consistent cell quality

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ISM

Performance Maintenance & Optimization

Optimal performance to yield top product quality:

- Professional installation
- Well documented traceable configuration and settings
- Proactive maintenance
- Protection or restoration of your equipment's reliability

Compliance Assurance & Certification

In full compliance with standards & regulations:

- Calibration and performance verification
- Audit-proof documentation of equipment qualification testing
- Recommendations for routine testing and maintenance schedules

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